

**KAMIYA BIOMEDICAL COMPANY**

# SOD Assay

**For the measurement of Superoxide Dismutase (SOD) Inhibition Activity**

**Cat. No. KT-019 (500 tests)**

**For Research Use Only. Not for Use in Diagnostic Procedures.**

**PRODUCT INFORMATION****SOD Assay**  
Cat. No. KT-019 (500 tests)**PRODUCT**

The **K-ASSAY**® SOD Assay is for the measurement of Superoxide Dismutase (SOD) inhibition activity.

**PRINCIPLE**

Superoxide Dismutase (SOD), which catalyzes the dismutation of the superoxide anion ( $O_2^-$ ) into hydrogen peroxide and molecular oxygen, is one of the most important antioxidative enzymes. In order to determine the SOD activity, several direct and indirect methods have been developed. Among these methods, an indirect method using nitroblue tetrazolium (NBT) is commonly used due to its convenience and ease of use. However, there are several disadvantages to the NBT method, such as poor water solubility of the formazan dye and the interaction with the reduced form of xanthine oxidase. The **K-ASSAY**® SOD Assay allows very convenient SOD assaying by utilizing **KAMIYA BIOMEDICAL COMPANY**'s highly water-soluble tetrazolium salt, WST-1 (2-(4-Iodophenyl)-3-(4-nitrophenyl)-5-(2,4-disulfophenyl)-2H-tetrazolium, monosodium salt) that produces a water-soluble formazan dye upon reduction with a superoxide anion. The rate of the reduction with  $O_2^-$  are linearly related to the xanthine oxidase (XO) activity, and is inhibited by SOD, as shown in Figure 1. Therefore, the  $IC_{50}$  (50% inhibition activity of SOD or SOD-like materials) can be determined by a colorimetric method.




Fig. 1 Principle of the determination of SOD activity using SOD Assay Kit.

## COMPONENTS

KT-019 (500 Tests)

- 1 x 5 mL WST Solution
- 1 x 100  $\mu$ L Enzyme Solution
- 1 x 100 mL Buffer Solution
- 1 x 50 mL Dilution Buffer

## MATERIALS OR EQUIPMENT REQUIRED BUT NOT PROVIDED

- Plate reader (450 nm filter).
- 96-well microplate.
- 2-20  $\mu$ L & 20-200  $\mu$ L pipettes and a multi-channel pipette.
- Incubator.

## PROTOCOLS

### Preparation of Solutions (for one 96-well plate)

- WST Working Solution  
Dilute 1 mL of WST Solution with 19 mL of Buffer Solution. WST working solution is stable for 2 months at 4°C.
- Enzyme Working Solution  
Centrifuge the Enzyme Solution tube for 5 sec. Mix by pipeting, and dilute 15  $\mu$ L of Enzyme Solution with 2.5 mL of Dilution Buffer. Enzyme working solution is stable at 4°C for 3 weeks.
- Sample Solution  
Dilute sample solution with Dilution buffer or saline to prepare sample solution as follows.  
Dilution ratio: 1, 1/5, 1/5<sup>2</sup>, 1/5<sup>3</sup>, 1/5<sup>4</sup>, 1/5<sup>5</sup>, 1/5<sup>6</sup>




Fig. 2 Preparation of sample solution

Table 1: Amount of each solution for sample, blank 1, 2 and 3

	sample	blank 1	blank 2*	blank 3
Sample Solution	20 $\mu$ L	-	20 $\mu$ L	-
ddH <sub>2</sub> O	-	20 $\mu$ L	-	20 $\mu$ L
WST Working Solution	200 $\mu$ L	200 $\mu$ L	200 $\mu$ L	200 $\mu$ L
Enzyme Working Solution	20 $\mu$ L	20 $\mu$ L	-	-
Dilution Buffer	-	-	20 $\mu$ L	20 $\mu$ L

**blank 1: coloring without inhibitor**

**blank 2: sample blank**

**blank 3: reagent blank**

**\*If the color of the sample solution is strong, measure blank 2 at each dilution of the sample.**

### General Protocol

1. Add 20  $\mu$ L of sample solution to each sample and blank 2 well, and add 20  $\mu$ L of ddH<sub>2</sub>O (double distilled water) to each blank 1 and blank 3 well.
2. Add 200  $\mu$ L of WST Working Solution to each well, and mix.
3. Add 20  $\mu$ L of Dilution Buffer to each blank 2 and blank 3 well.
4. Add 20  $\mu$ L of Enzyme Working Solution to each sample and blank 1 well, and then mix thoroughly\*\*.
5. Incubate the plate at 37°C for 20 min.
6. Read the absorbance at 450 nm using a microplate reader.
7. Calculate the SOD activity (inhibition rate %) using the following equation.

$$\text{SOD activity (inhibition rate \%)} = \{[(A_{\text{blank 1}} - A_{\text{blank 3}}) - (A_{\text{sample}} - A_{\text{blank 2}})] / (A_{\text{blank 1}} - A_{\text{blank 3}})\} \times 100$$

\*\* Since superoxide will be released immediately after the addition of Enzyme Working Solution to a well, use a multi-channel pipette to avoid the reaction time lag of each well.

Inhibition activity can also be determined by a kinetic method. Please determine an incubation time range that has a linearity of the slope before the assay. A good linearity should be observed up to 20 min. For the calculation, use the following equation:

$$\text{SOD activity (inhibition rate \%)} = \{[(S1 - S3) - (SS - S2)] / (S1 - S3)\} \times 100$$

S1: slope of blank1, S2: slope of blank2, S3: slope of blank3, SS: slope of sample

	1	2	3	4	5	6	7	8	9	10	11	12
A	Sample 1	1										
B	Sample 1	1/5										
C	Sample 1	1/5 <sup>2</sup>										
D	Sample 1	1/5 <sup>3</sup>	Sample 2			Sample 3			Sample 4			
E	Sample 1	1/5 <sup>4</sup>										
F	Sample 1	1/5 <sup>5</sup>										
G	Sample 1	1/5 <sup>6</sup>										
H	blank 1	blank 2			blank 3							

Fig 3. Example of sample and blank arrangement on a 96-well plate




Fig 4. Inhibition curve of Cu,Zn-SOD (SIGMA Cat#:S2515)

One unit of SOD is defined as the amount of the enzyme in 20  $\mu\text{L}$  of sample solution that inhibits the reduction reaction of WST-1 with superoxide anion by 50%.

### Determination of SOD Activity

1. Read the dilution ratio at 50% inhibition ( $\text{IC}_{50}$ ) from inhibition curve.
2. Multiply the dilution ratio at  $\text{IC}_{50}$  and at the sample preparation to obtain the SOD activity.

Example: Determination of SOD activity in erythrocytes

1. Read the dilution ratio at 50% inhibition ( $\text{IC}_{50}$ ) from the inhibition curve (Fig.5). The dilution ratio at  $\text{IC}_{50}$  is 1/1.8.
2. The SOD activity before dilution is "1.8 U/20  $\mu\text{L}$ " from the definition above.
3. The SOD activity in 1 mL of sample is  $1.8 / 0.02 = 90.0 \text{ U/mL}$ .
4. Multiply the SOD activity calculated above by the necessary dilution ratio for sample preparation. In the case blood is diluted to 1/108 during the sample preparation, the SOD activity in blood is calculated as below.  
 $108 \times 90.0 = 9,720 \text{ U/mL}$  of blood.




Fig.5 Inhibition curve of erythrocyte

### Determination of Mn-SOD activity

Mn-SOD activity can be measured by adding potassium cyanide (final concentration: 1 mmol/L) or diethyldithiocarbamate (final concentration: 1 mmol/L) to the sample solution. These reagents inactivate Cu, Zn-SOD and extracellular-SOD activities.

### INTERFERENCE

Below shows compatible concentration of possible interfering materials. If sample contains these materials, please dilute the sample to be below their compatible concentration. Since 2-mercaptoethanol and dithiothreitol cause a significant increase of the O.D. value, please remove them when sample contains these materials.

Detergents: SDS (0.05%), Tween 20 (0.5%), NP-40 (0.5%)

Solvents: Ethanol (25%), DMSO (5%)

Reducing agents: Glutathione reduced form (1.25 mmol/L), Ascorbic acid (0.1 mmol/L)

Others: EDTA (2 mM), BSA (1%w/v)

### PRECAUTION

1. For dilution of sample, use Dilution buffer or saline.
2. The Enzyme solution is separated into two layers. Therefore, omitting the pipetting process will result in inaccurate experimental results.
3. For an accurate measurement, the use of multiple wells per sample is recommended.
4. Since superoxide will be released immediately after the addition of Enzyme working solution to a well, use a multi-channel pipette to avoid the reaction time lag of each well.

### STORAGE

Please store at 4°C. Protect the WST Solution and WST Working Solution from light.

**FOR RESEARCH USE ONLY  
NOT FOR USE IN DIAGNOSTIC PROCEDURES**

## **KAMIYA BIOMEDICAL COMPANY**

12779 Gateway Drive, Seattle, WA 98168  
Tel: (206) 575-8068 Fax: (206) 575-8094  
Email: LifeScience@k-assay.com  
www.k-assay.com