


Animal CRP ELISA Kits

C-Reactive protein (CRP) is produced in the liver and levels of CRP rise in response to inflammation in the body. This acute-phase protein binds to lysophosphatidylcholine expressed on dead or dying cells in order to activate the complement immune system.



NEW KITS:

Product	Cat. No.	Sample Type	Assay Range
Bovine CRP ELISA	KT-784	Serum, plasma, and milk	0.98 - 62.5 ng/mL
Goat CRP ELISA	KT-785	Serum and plasma	3.91 - 125 ng/mL
Sheep CRP ELISA	KT-786	Serum and plasma	7.81 - 250 ng/mL

Classic Kits:

Product	Cat. No.	Sample Type	Assay Range
Dog CRP ELISA*	KT-093	Serum and plasma	3.125 - 200 ng/mL
Mouse CRP ELISA*	KT-095	Serum and plasma	0.39 - 25 ng/mL
Rabbit CRP ELISA*	KT-097	Serum and plasma	6.25 - 200 ng/mL
Rat CRP ELISA*	KT-099	Serum and plasma	6.25 - 200 ng/mL
Pig CRP ELISA*	KT-184	Serum and plasma	6.25 - 100 ng/mL
Monkey CRP ELISA	KT-425	Serum	1.17 - 75 ng/mL
Cat CRP ELISA	KT-486	Serum and plasma	3.125 - 200 ng/mL
Horse CRP ELISA	KT-487	Serum and plasma	10 - 320 ng/mL

* Includes Positive Control!



For *In Vitro* Research Use Only. Not for use in diagnostic procedures in the U.S.
 For additional information please visit www.K-ASSAY.com

KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive, Seattle, WA 98168
 Tel: (206) 575-8068 Fax: (206) 575-8094
 Email: LifeScience@k-assay.com
www.k-assay.com