




KL-6 (Krebs von den Lungen-6) Assay Reagent for Chemistry Analyzers

Fully automated immunoturbidimetric assay for the quantitative measurement of KL-6 in serum and plasma samples on chemistry analyzers.

KL-6 is a mucinous high-molecular weight glycoprotein. Research studies have reported that KL-6 levels are significantly higher in individuals with active pneumonitis, and that KL-6 is an important biomarker for interstitial lung diseases.

Features:

- **Uses a specific monoclonal antibody to KL-6**
- **Liquid, ready to use reagents**
- **Can be used with serum or plasma samples**
- **Assay Range: 70 - 10,000 U/mL**
- **Applications available for most chemistry analyzers.
Contact us for more information.**



Ordering Information:

| | |
|-----------------|---|
| KAI-250 | KL-6 Reagent (R1: 1 x 15 mL, R2: 1 x 5 mL) |
| KAI-251C | KL-6 Calibrator (5 x 1 mL, 5 concentrations) |
| K252C-2M | KL-6 Control (2 levels: 1 x 2 mL each level) |

For research use only. Not for use in diagnostic procedure in the U.S.

KAMIYA BIOMEDICAL COMPANY

12779 Gateway Drive, Seattle, WA 98168

800-526-4925 • 206-575-8068 • FAX: 206-575-8094

www.k-assay.com