



## PRODUCT DATA SHEET

**Product:** Anti-Aggrecan mAb, clone 5B11

**Cat. No.:** MC-1004 (100 µL)

**Specificity:**

Reacts with the conformational epitope of aggrecan.

**Species Reactivity:**

Human. Other species not tested.

**Ig Isotype:**

Mouse IgG1

**Immunogen:**

Human articular cartilage aggrecan

**Format:**

100 µL of ascites.

**Storage:**

Store at 4°C for short term and at -20°C for long term. Aliquot to avoid freeze/thaw cycles.

**Applications:**

- ELISA
- Immunohistochemistry
- Immunoprecipitation
- Western blot

The optimal dilution for a specific application should be determined by the researcher.

**Limitations:**

For *in vitro* research use only. Not for use in diagnostics or in humans.

**Warranty:**

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.