



PRODUCT DATA SHEET

Product: Anti-human-CD40Ligand mAb, clone MK13A4

Cat. No.: MC-277 (100 µg)

Specificity:

This antibody recognizes CD40Ligand.

Species Reactivity:

Human, Others not tested.

Ig Isotype:

Mouse IgG1

Format:

100 µg/100 µL of purified monoclonal antibody in PBS.

Storage and Stability:

Store at 4°C. After opening store at -20°C in aliquots. Avoid repeat freeze/thaws.

Applications and Suggested Dilutions:

- Functional assays: This antibody was found to block the binding of CD40 Ligand to its receptor CD40.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.