



PRODUCT DATA SHEET

Product: Anti-mouse CD40Ligand, Clone MR1

Cat. No: MC-278 (100 µg)

Specificity:

This antibody reacts with CD40Ligand (gp39).

Species Reactivity:

Mouse, Others not tested.

Ig Isotype: Hamster IgG1

Format:

100 µg/100 µl of purified antibody in PBS.
Protein A purification.

Storage and Stability:

Unopened store at 2-8°C. After opening aliquot
and store at -20°C.

Applications and Suggested Dilutions:

■ *in vivo* studies

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.