



PRODUCT DATA SHEET

Product: Anti-DR6 mAb, clone DR-6-04-EC

Cat. No.: MC-345 (100 µg)

Synonyms:

Anti-Death Receptor-6 (human) MAb (DR-6-04-EC), anti-TNFRSF 21 (human) MAb (DR-6-04-EC)

Specificity:

Recognizes human DR6.

Species Reactivity:

Human. Others not tested.

Ig Isotype:

Mouse IgG₁

Immunogen:

Recombinant human DR6 (death receptor-6) extracellular domain (aa 42-335).

Format:

100 µg of a 1 mg/mL of protein A-affinity purified antibody in PBS, pH 7.4, containing 15 mM sodium azide. Purity >95% by SDS-PAGE.

Storage:

Store at 4°C, short term. Store at -20°C long term. Aliquot to avoid freeze/thaw cycles.

Applications:

- Flow cytometry
- Immunocytochemistry
- Immunoprecipitation

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.