



PRODUCT DATA SHEET

Product: Anti-TIRAP, clone Sebi-1

Cat. No.: MC-394 (100 µg)

Synonyms:

Toll-interleukin 1 Receptor Domain-containing Adaptor Protein, Mal, MyD88 Adapter-like Protein.

Specificity:

Recognizes human TIRAP.

Species Reactivity:

Human. Other species not tested.

Ig Isotype:

Rat IgG_{2a}

Immunogen:

Recombinant human TIRAP (aa 1-215 of TIRAP was fused at the N-terminus to GST).

Format:

100 µg of antibody purified from concentrated hybridoma tissue culture supernatant at 1 mg/mL in PBS buffer with 0.02% sodium azide. ≥95% purity by SDS-PAGE.

Storage:

Store at 4°C for short term and at -20°C long term. Aliquot to avoid freeze/thaw cycles.

Applications and Suggested Dilutions:

- Western blot: Use at a 1:500 dilution.
- Immunoprecipitation: Use at a 1:200 dilution.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.