

PRODUCT DATA SHEET

Product: Anti-Glutathione mAb, clone D8

Cat. No.: MC-750 (100 µg)

Specificity:

Recognizes glutathione.

Host:

Mouse

Ig Isotype:

IgG2a

Immunogen:

KLH conjugated glutathione

Format:

100 µL of 1 mg/mL anti-glutathione antibody in 1X PBS (pH 7.2), 0.01% sodium azide. Protein A purified.

Storage and Stability:

4°C short term, -20°C long term. Aliquot to avoid cycles of freeze/thaw.

Applications and Suggested Dilutions:

- ELISA: Use at a 1:1,000 dilution.
- Immunohistochemistry: Use at a 1:100 dilution.
- Immunoprecipitation
- Western blotting: (non-reducing conditions), Use at a 1:1,000 dilution.

The optimal dilution for a specific application should be determined by the researcher.



Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.