



PRODUCT DATA SHEET

Product: Anti-Rabbit IgG Negative Control pAb

Cat. No.: PC-124 (500 µg)

Background:

This IgG was purified from serum collected from a rabbit prior to immunization. It does not react with any antigen of human or animal cells. In order to serve as a true negative control, this IgG should be used under the identical conditions and at the same concentration as KAMIYA BIOMEDICAL COMPANY's other purified rabbit polyclonal IgG antibodies.

Format:

500 µg of 1 mg/mL polyclonal antibody in 10 mM PBS, pH 7.4, with protein stabilizer and 0.09% sodium azide. Purified from rabbit serum by Protein A chromatography.

Storage:

Store at 4 °C.

Applications:

Negative control for:

- Flow cytometry
- ELISA
- Gel supershift
- Immunohistochemistry (Frozen and formalin / paraffin)
- Immunofluorescence
- Immunoprecipitation
- Western blot

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.