

PRODUCT DATA SHEET

Product: Anti-Hemagglutinin (HA) pAb

Cat. No.: PC-259 (0.1 mg)

Procedure:

Rabbits were immunized with HA peptide conjugated to KLH. After multiple immunizations in Freund's adjuvant, serum was collected. Antibodies were immunoaffinity-purified using the peptide conjugated to a solid-phase support.

Host:

Rabbit

Format:

Liquid. 0.1 mL of affinity-purified antibody at 1.0 mg/mL in PBS pH 7.2 and 0.1% sodium azide.


Storage:

Store at 4°C.

Applications and Suggested Dilutions:

- ELISA: 1.0 – 4.0 µg/mL (capture); 0.05 – 0.2 µg/mL (detection)
- Immunofluorescence: 1.0 – 4.0 µg/mL
- Western blot: 0.2 – 1.0 µg/mL

The optimal dilution for a specific application should be determined by the researcher.



Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.