



PRODUCT DATA SHEET

Product: Anti-Hexahistidine pAb

Cat. No.: PC-276 (1.0 mg)

Procedure:

Rabbits were immunized with HHHHHH peptide conjugated to KLH. After multiple immunizations in Freund's adjuvant, serum was collected. Antibodies were immunoaffinity-purified using the peptide immobilized on a solid support.

Host:

Rabbit

Format:

Liquid. 1.0 mL of affinity-purified polyclonal antibody in PBS, pH 7.2 and 0.1% sodium azide.

Storage:

Store at 4°C.

Applications and Suggested Dilutions:

- ELISA: Use at a 1:200 – 1:500 dilution for coating. Use at a 1:2,000 – 1:30,000 dilution for detection.
- Western blot: Use at a 1:1,000 – 1:3,000 dilution.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.