



PRODUCT DATA SHEET

Product: Anti-C-Myc (EQKLISEEDL), HRP-labeled, Polyclonal

Cat. No: PC-311 (0.1 mg)

Procedure:

Rabbits were immunized with the C-Myc peptide conjugated with KLH. After multiple immunizations in Freund's adjuvant, serum was collected. Antibodies were immunoaffinity-purified using the peptide immobilized on a solid phase and conjugated with horseradish peroxidase (HRP).

Host: Rabbit

Format:

Liquid. HRP-conjugated antibody in 50 mM HEPES pH 7.0, 0.1 M sodium chloride with 0.2% BSA and 0.01% thimerosal.

Storage and Stability:

Store at +4°C.

Applications and Suggested Dilutions:

- ELISA: coating 1:5,000 – 1:20,000
- Western blot: Use at 1:1,000-1:5,000.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.