



PRODUCT DATA SHEET

Product: Anti-SET, polyclonal

Cat. No.: PC-381 (50 µg)

Synonyms:

HLA-DR-associated Protein 2, PHAPII, PP2A Inhibitor Protein 2 α , I 2α ^{PP2A}, IGADD, Inhibitor of Granzyme-activated DNase.

Specificity:

Recognizes human SET. Detects endogenous protein by IP and WB.

Species Reactivity:

Human. Other species not tested.

Immunogen:

Synthetic peptide corresponding to aa 3-16 (PKRQSPLPPQKKKP) of human SET.

Host:

Rabbit

Format:

50 µg of protein A purified antibody at 1 mg/mL in PBS and 0.02% sodium azide.

Storage:

Store at 4°C short term, -20°C long term. Aliquot to avoid repeat freeze/thaw cycles.

Applications and Suggested Dilutions:

- Immunocytochemistry
- Immunoprecipitation
- Flow cytometry: Use at a 1:100-1:200 dilution.
- Western blots: Use at a 1:1,000 dilution.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.