



## PRODUCT DATA SHEET

**Product:** Anti-Adiponectin Receptor 2 pAb

**Cat. No.:** PC-450 (100 µL)

**Synonym:**  
Anti-AdipoR2

**Specificity:**  
Recognizes mouse Adiponectin Receptor 2.

**Species Reactivity:**  
Mouse. Other species not tested.

**Immunogen:**  
Synthetic peptide corresponding to aa 2 – 36 (N<sup>2</sup>EPAKHRLGCTRTPEPDIRLRKGHQLDDTRG SNND<sup>36</sup>) of mouse Adiponectin Receptor 2.

**Host:**  
Rabbit

**Format:**  
100 µL of neat liquid mouse serum containing 0.02% sodium azide.

**Storage:**  
Store at 4°C short term, -20°C long term. Aliquote to avoid freeze/thaw cycles.

**Application:**  
■ Immunoprecipitation: Use at a 1:50 dilution.  
The optimal dilution for a specific application should be determined by the researcher.

**Limitations:**  
For *in vitro* research use only. Not for use in diagnostics or in humans.

**Warranty:**  
No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.