



PRODUCT DATA SHEET

Product: Anti-CTRP7, polyclonal

Cat. No.: PC-485 (50 µg)

Synonyms:

C1q Tumor Necrosis Factor- α -related Protein 7;
C1qTNF7

Specificity:

Recognizes mouse CTRP7.

Species Reactivity:

Mouse. Other species not tested.

Host:

Rabbit

Immunogen:

Recombinant mouse CTRP7 (aa 18-290).

Format:

50 µg of epitope-affinity purified antibody at 1 mg/mL in PBS with 0.02% sodium azide.

Storage:

Store at 4°C short term, -20°C long term. Aliquot to avoid repeat freeze/thaw cycles.

Applications and Suggested Dilutions:

- Western blot: Use at a 1:1,000 dilution
- Immunoprecipitation: Use at a 1:200 dilution

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.