



PRODUCT DATA SHEET

Product: Anti-TWEAK pAb

Cat. No.: PC-497 (50 µg)

Synonyms:

Apo-3L; TNFSF 12

Specificity:

Recognizes human and mouse TWEAK.

Species Reactivity:

Human and Mouse. Other species not tested.

Host:

Rabbit

Immunogen:

Recombinant human TWEAK. The extracellular domain of human TWEAK (aa 106-249).

Format:

50 µg of epitope-affinity purified polyclonal antibody at 1 mg/mL in PBS containing 0.02% sodium azide.

Storage:

Store at 4°C short term, -20°C long term. Aliquot to avoid repeat freeze/thaw cycles.

Applications and Suggested Dilutions:

- Western blot: Use at a 1:1,000 dilution.
- Immunoprecipitation: Use at a 1:100 dilution.
- Flow cytometry: Use at a 1:1,000 dilution.

The optimal dilution for a specific application should be determined by the researcher.

Limitations:

For *in vitro* research use only. Not for use in diagnostics or in humans.

Warranty:

No warranties, expressed or implied, are made regarding the use of this product. KAMIYA BIOMEDICAL COMPANY is not liable for any damage, personal injury, or economic loss caused by this product.